


scimaX

● Identify with Confidence

Taking science to the max

Remarkable advances have been made in life sciences research using the advanced capabilities of the versatile solarix platform. Now imagine if this enabling technology was made smaller and maintenance-free when it came to liquid cryogenics, so that it could fit in a standard laboratory.

scimaX allows you to dramatically improve your productivity by operating around the clock by doing your ESI experiments during the day and acquiring MALDI imaging data when you go home. Take on projects that demand exceptional scientific insight and breakthrough with scimaX.




A maximum resolving power of > 20,000,000 without the need for operational cryogen fills, ever!

Solve difficult problems in science and capture the "high hanging fruit." When you want the ultimate extreme performance, with advanced versatility, the only choice is scimaX.

- Featuring an integrated 7T conduction cooled magnet, scimaX is small and fits in a standard lab
- Liquid helium filling and quench line not required
- Includes MALDI and ESI with no switching required
- Continuous Accumulation of Selected Ions (CASI) for enhanced sensitivity
- Choice of 8 fragmentation techniques


Footprint: scimaX and solarix




solarix

eXtreme Resolution MALDI Imaging powers new discoveries

scimaX is the ultimate MALDI imaging system for analyzing small to medium molecules, m/z 100-1,500. Its unrivaled eXtreme Resolution capability and sub-ppm mass accuracy, over a wide mass range, can differentiate images that are only mDa apart and are prerequisite for Isotopic Fine Structure (IFS) analysis and molecular formula confirmation.




Realize the full benefit of eXtreme Resolution imaging. In the 3-color image, above, the distribution of 3 ions that have a nominal m/z of 848 and differing by only 16 mDa are spatially differentiated.



scimaX Imaging of a rat brain with standard 2xR detection provides fast imaging speed while retaining eXtreme mass resolution.

Mass accuracy and Isotopic Fine Structure (IFS) are indispensable when determining a unique molecular formula for an image and, just as importantly, to differentiate isobaric interfering peaks that bias distribution measurements.






scimaX - powered by MRMS technology

Isotope	mass (Da)	abundance(%)
¹ H	1.007825	99.9885
² H	2.014102	0.0115
¹² C	12.000000	98.9300
¹³ C	13.003355	1.0700
¹⁴ N	14.003074	99.6360
¹⁵ N	15.000109	0.3640
¹⁶ O	15.994915	99.7570
¹⁷ O	16.999132	0.0380
¹⁸ O	17.999161	0.2050

Isotope	mass (Da)	abundance(%)
³² S	31.972071	94.9900
³³ S	32.971459	0.7500
³⁴ S	33.967867	4.2500
³⁶ S	35.967081	0.0100
³⁵ Cl	34.968853	75.7600
³⁷ Cl	36.965903	24.2200
³⁹ K	38.963707	93.2581
⁴⁰ K	39.963998	0.0117
⁴¹ K	40.961826	6.7302

$m(e^-) = 0.000549$ Da






MRMS aXelerate – high throughput with unmatched specificity for Phenomics research

- Accelerate sample throughput enabling large cohort and longitudinal studies in phenomics research (> 200 samples/day)
- Simultaneous analysis of known and unknown metabolites
- Access compounds not readily detectable by LC-MS analysis

Whether Metabolomics, Phenomics or any other complex sample analysis, large scale sample evaluation is now possible as MRMS aXelerate utilizes the eXtreme resolving power of the scimaX MRMS and **MetaboScape 4.0** to enable a powerful LC free solution.

eXtreme Resolution (XR) allows for direct sample analysis enabling true high sample throughput complementary to established NMR based solutions. From the largest unknown to the smallest, MRMS **aXelerate** incorporates a combination of ultra-high mass accuracy, True Isotopic Pattern, and Isotopic Fine Structure, allowing confident assignments of molecular formulae at any level.

Flow Injection Analysis (FIA) or MALDI based workflows provide access to compounds not readily detectable by LC-MS and provide high information content across a wide dynamic range.



1
Highest abundant compound annotated as: C₆H₁₂O₆Na
mass accuracy:
0.09 ppm
RP: 1,500,000

2
Expanded A+1 region of C₆H₁₂O₆Na reveals Isotopic Fine Structure for this compound. Such hexose sugars are not well retained on reversed phase LC-MS and difficult to detect.

3
Low intensity ion with 1.1 × 10⁶ vs. highest abundant peak with 3.7 × 10⁹ intensity demonstrate > 3 order dynamic range.

4
Mass region zoom: 267.00 – 267.40 m/z. This 0.4 Da mass region exemplifies the richness of information in FIA-MRMS data: 6 species annotated.

scimaX

scimaX is easy to site and maintain

Fits in a standard lab, no quench line needed, no filling of liquid helium needed

scimaX is an integrated and versatile instrument

ESI and MALDI sources are standard and compatible with various API sources (APPI, APCI, GC-APCI) and includes several ion activation techniques (CID, (n)ETD, (n)ECD, EID, SORI-CID, MALDI-MSD)

scimaX brings "high hanging fruit" within reach ...

... for various applications like MALDI Imaging, Phenomics, Petroleomics, and more

scimaX enables identification with confidence

Routine sub-ppm mass accuracy and Isotopic Fine Structure (IFS) capability enable high-confidence molecular formula assignments for known and unknown targets



**Professor Evgeny Nikolaev,
ParaCell Inventor,
Skoltech, Moscow**

"With scimaX, Bruker continues innovation of MRMS technology to reduce effort while increasing scientific output. scimaX' core detection technology, the 2xR ParaCell, provides uncommon broadband ion stability and mitigates ion cloud coalescence resulting in mass resolution orders of magnitude above other detection schemes. This enables extreme resolving power over a broad mass range needed to perform isotopic fine structure analysis of complex mixtures."

US7315020B2, GB2402261B, US7064321B2, DE10213652B4, GB2390937B, US6803569B2, DE102009050039B4, US8704173B2, US8859953B2, US8766174B1, DE102014226498B4, US9620349B2, EP2858090B1, US9355830B2, US9111735B1

For research use only. Not for use in diagnostic procedures.

● **Bruker Daltonics GmbH & Co. KG**

Bremen · Germany
Phone +49 (0)421-2205-0

Bruker Scientific LLC

Billerica, MA · USA
Phone +1 (978) 663-3660



**Scan the QR-Code
for more Details**

ms.sales.bdal@bruker.com - www.bruker.com