

MACROMOLECULAR CRYSTALLOGRAPHY


D8 VENTURE D8 QUEST

High-end Solutions for
X-ray Crystallography

Innovation with Integrity




IµS DIAMOND II: Rotating anode performance




SCOUT: Automated cryo-cooled protein crystal handling



Low temperature attachment




ISX Stage: Automated *in-situ* crystal screening



Stable and easy-to-align beam path



Wide doors for easy access



METALJET D2 PLUS: the most intense home source



PHOTON III photon-counting detector with large active area





X-ray Crystallography drives Structural Biology forward

The field of structural biology continues to advance dramatically with researchers tackling ever more complex and challenging macromolecular structures. X-ray crystallography, the engine of this ongoing quest, is accounting for over 75% of the structures deposited in the Protein Data Bank since 2020.

The best in-house X-ray systems provide the fastest and most economical route to protein structures for increasingly challenging problems, such as large protein complexes and tiny crystals.


Great research deserves great tools.

As structural biology advances, the need for fast, economical structure determination is becoming more acute. Modern laboratories, whether in academic research or commercial drug-discovery, need better crystallographic instruments in-house – more powerful, faster and much easier to use.


At Bruker, our driving passion is to provide scientists with highest-quality tools for crystallography. With this goal in mind, our research and development teams are continuously pushing the frontiers of innovation, helping you get better structures faster with greatly increased efficiency. These efforts resulted in a range of highly advanced home-lab instruments:

Third-generation Macromolecular Crystallography Solutions – D8 QUEST and D8 VENTURE.


Metalloproteins ¹⁾




Molecular motors ⁴⁾




High-throughput crystallography ⁷⁾




High-resolution structures ¹⁰⁾




GPCR crystallography ¹³⁾




Microcrystals ²⁾




Structural proteins and scaffolds ⁵⁾




Protein complexes ⁸⁾




Fragment-based lead generation ¹¹⁾




Protein-protein interactions ¹⁴⁾




SAD phasing ³⁾




Multi-protein complexes ⁶⁾




Membrane proteins ⁹⁾



Molecular replacement ¹²⁾



Viral protein-RNA structures ¹⁵⁾




1) Kostan, J; Sjöblom, B; Maixner, F; Mlynek, G; Furtmüller, P G; Obinger, C; Wagner, M; Daims, H; Djinić-Carugo: *J. Struct. Biol.* 2010, **172**, 331.
 2) Benning, M: Bruker application laboratory Madison, WI, USA.
 3) Beck, T; Krasauskas, A; Gruene, T; Sheldrick, G M: *Acta Cryst.* 2008, **D64**, 1179.
 4) Hernandez-Valladares, M; Kim, T; Kannan, B; Tung, A; Aguda, A H; Larsson, M; Cooper, J A; Robinson R C: *Nat. Struct. Mol. Biol.* 2010, **17**, 497.

5) Cleverley, R M; Rutter, Z J; Rismondo, J et al.: *Nat Commun.* 2019, **10**, 261.
 6) Wang, F; Mei, Z; Qi, Y; Yan, C; Hu, Q; Wang, J; Shi, Y: *Nature* 2011, **471**, 331. *)
 7) Benning, M; Biadene, M: Bruker application laboratory Madison, WI, USA.
 8) Bochtler, M; Ditzel, L; Groll, M; Hartmann, C; Huber, R: *Annu. Rev. Biophys. Biomol. Struct.* 1999, **28**, 295.
 9) Prince, S M; Achtman, M; Derrick, J P: *PNAS* 2002, **99**, 3417.

10) Benning, M: 2009, Lab Report SC-XRD 42, DOC-L86-EXS042.
 11) Hoarau, M; Vanichtanukul, J; Srimongkolpithak, N; Vitsupakorn, D; Yuthavong, Y; Kamchonwongpaisan, S: *J. Enzyme Inhib. Med. Chem.* 2021, **36**, 198.
 12) Bolanos-Garcia, V M, Chirgadze, D Y; Blundell, T L: 2011 BUBR1 kinase, personal communication.
 13) Doré, A. S. (Heptares Therapeutics Ltd.); Freisz, S. (Bruker AXS): Data Collection on D8 VENTURE with METALJET; 2016.

14) Makraki, E; Darby, J F; Carneiro, M G; Firth J D; Heyam, A; Ab, E; O'Brien, P; Siegal, G; Hubbard, R E: *Biochem.*, 2020, **477**, 4383.
 15) Krafcikova, P; Silhan, J; Nencka, R et al: *Nat Commun.*, 2020, **11**, 3717.

*) Plot created using data from the Protein Data Bank: Berman, H M; Westbrook, J; Feng, Z; Gilliland, G; Bhat, T N; Weissig, H; Shindyalov, I N; Bourne P E: *Nucleic Acids Res.* 2000, **28**, 235.



IµS DIAMOND II X-ray source with HELIOS optics

Brighten up your home lab with IµS DIAMOND II or METALJET D2 PLUS


High flux, stability, small beam size, and low divergence are essential for collecting high-quality data from challenging crystals. These requirements have driven developments in synchrotron design. Now, Bruker has applied these important principles to home-lab source technology. Our sources deliver beam intensities comparable to those of typical bending-magnet beamlines. The modern METALJET D2 PLUS crosses the next boundary, matching diffraction limits previously only seen at third-generation synchrotrons. These solutions combine high intensity and superior beam profiles with unmatched beam stability. Instant instrument access makes home-lab crystallography with Bruker solutions more productive than ever before.

Copper is king.

Copper radiation (1.54 Å) is the most popular wavelength for macromolecular crystallography in the home lab. X-rays at this wavelength, provided by our IµS DIAMOND II, are diffracted very strongly by the light atoms in biological molecules. At Cu radiation's wavelength, sulfur and other common scatterers also provide the anomalous signal for in-house SAD phasing.

Gallium: the better copper.

Home Sources vs. Bending Magnet Beamlines



The gallium-K α emission line (1.34 Å) is near that of Cu for maintaining the proven benefits of Cu-K α radiation. Yet, Gallium's shorter wavelength adds a number of advantages, such as reduced radiation damage, increased data collection efficiency, reduced X-ray absorption by the sample, and reduced scattering by both mother liquor and air. The METALJET's ultra-intense Ga-K α beam takes advantage of all these benefits.

Your personal beamline.

In the past, finding the right X-ray source for your lab demanded a choice between the highest performance and the lowest maintenance. Now you can have it all: Bruker's new X-ray sources are designed to combine superb performance with unprecedented stability, maximum uptime, and easy maintenance. All our compact X-ray sources are directly mounted on the goniometer. This ensures easy and stable downstream alignment of the source to the goniometer – greatly improving overall system precision and minimizing service efforts. All generators are fully contained within the cabinet, resulting in very compact systems with minimal footprints.

IPUS DIAMOND II: sealed tube convenience, rotating anode performance.

With over 1,700 systems installed in the 20 years since its introduction, the air-cooled IPUS has developed a formidable reputation as the most reliable low-power microfocus X-ray source available, generating an impressively intense X-ray beam.

All our IPUS sources are developed completely and specifically for crystallography which results in higher intensity than any other microfocus source. However, there is more we can do: diamond hybrid anode technology, with an eight times better heat conductivity than copper, allows a further increase in the power load, leading to the IPUS DIAMOND, a source outperforming modern microfocus rotating anodes.

The second generation, the IPUS DIAMOND II, pushes the limits in microfocus source technology for biological crystallography even further. The IPUS DIAMOND II features isotopically pure diamond, which has the highest thermal conductivity of any known material. A new high-brightness cathode produces a more homogenous, higher current density electron beam. Both improvements allow a power load increased once more, leading to ultimate X-ray intensities and making conventional microfocus rotating anode generators obsolete. Even better, the IPUS DIAMOND II has essentially no maintenance costs or downtime. We support every single unit with a three-year warranty.

IPUS DIAMOND II – simply brilliant.

METALJET D2 PLUS with our best HELIOS optics

The X-ray revolution without rotation.

The METALJET D2 PLUS uses a high-speed jet of liquid metal that accepts a much higher power load than the solid target of rotating anodes. This breakthrough in X-ray source technology delivers an X-ray beam an order of magnitude brighter than that achieved with traditional rotating anode sources.

With the metal jet continuously supplying fresh target material, the source intensity remains constant over time, in contrast to rotating anodes that degrade in intensity quickly due to anode roughening. Like the I μ S DIAMOND II, the METALJET D2 PLUS is operated with single-phase power and also requires no external cooling, minimizing installation requirements.

HELIOS MX: the brightest optics under the sun.

Our HELIOS multilayer optic delivers up to three times the intensity of conventional multilayer optics. Improved deposition technology yields higher precision and greatly enhances the optic's reflectivity.

For the METALJET D2 PLUS, we specifically developed a unique low-figure-error optic. This patented optic is the only X-ray mirror that fully preserves the inherent brilliance of the METALJET source. Conversely, a METALJET combined with a conventional, relatively high-figure-error multilayer optic typically leads to a loss of brilliance of up to 70%.

Collimators for beam optimization.

The highly-monochromatic X-ray beam provided by the HELIOS optics can be easily and reproducibly optimized to match the characteristics of each sample. A set of collimators control beam size and divergence and is an integral part of each D8 QUEST and D8 VENTURE system.



Liquid metal METALJET D2 PLUS with HELIOS optics

Our X-rays are “Green”



Future-Proof Solutions

With the D8 QUEST and D8 VENTURE, we are following design principles that protect the environment. Our PHOTON III detector and advanced μ S DIAMOND II and METALJET D2 PLUS sources all feature very low energy consumption and do not need water cooling.

This significantly reduces carbon emissions and water consumption, improves reliability, and reduces the cost of ownership.

Saves 26 t of CO₂ per year *



No Water Supply

- Air-cooled PHOTON III
- Air-cooled μ S DIAMOND II microfocus source
- METALJET D2 PLUS, no external cooling water required

Saves 1,700 m³ cooling water per year *




Single-Phase Power

- Ultra-low power consumption of all PHOTON III detectors
- μ S DIAMOND II uses 99% less energy than conventional rotating anodes
- METALJET D2 PLUS delivers beam intensities more than three times higher than microfocus rotating anodes, using only a fraction of a rotating anode's power consumption
- Single-phase power and standard circuit breaker for ease of installation

Saves 43,000 kWh electrical energy per year *

*) 24 hours, 300 days with μ S DIAMOND II and PHOTON III




An entire tube production line – for the best μ S sources


All μ S tubes are designed, developed and manufactured exclusively by Incoatec. Indeed, Incoatec is the only company that makes microfocus tubes optimized for X-ray diffraction. These optimized tubes offer higher performance and higher quality, resulting in longer tube lifetimes.

Source	Typical tube lifetime	Relative intensity
Bruker μS DIAMOND II	5 years	1.1
Microfocus rotating anode	1 year	1.0 ^{*)}

^{*)} with fresh anode only

In contrast to rotating anodes, the heat load on the μ S DIAMOND II surface is constant over time, making the output of the μ S DIAMOND II much more stable. The μ S DIAMOND II source also outperforms microfocus rotating anodes without the headache of high running costs and without routine maintenance.





IµS DIAMOND II for Cu radiation


- The average intensity is 20% higher than that of a microfocus rotating anode source. The best optics put all the X-rays on the sample for up to 10 times lower scattered X-ray background.
- Unique 99% uptime guarantee.
- 10 times better stability than rotating anodes – for the best data quality.



METALJET D2 PLUS with HELIOS for Ga radiation

- The METALJET D2 PLUS uses a liquid target consisting of a Gallium rich alloy.
- Compared to Cu radiation, Ga-K α radiation allows higher multiplicity data sets in less time.
- Fresh, self-healing target for high power load with long-term constant beam intensity.


IµS DIAMOND II hybrid anode



Higher intensities are achieved in the IµS DIAMOND II by optimizing all the critical tube parts, such as the cathode, electron optics and take-off angles.

In addition, the diamond hybrid anode used in the IµS DIAMOND II features an isotopically pure diamond substrate that is coated with a layer of the target material and thus takes advantage of the high thermal conductivity of diamond to achieve even higher intensities.

The high-brightness cathode provides the most homogeneous electron-beam with the highest intensity of any microfocus tube.



The PHOTON III is available in three different sizes, perfectly matching the requirements of your application.

PHOTON III: Large active-area, photon-counting detector – no photon is left behind

The best crystal structures for publication require large active-area X-ray detectors with high signal-to-noise ratios for accurate weak and strong reflection intensities.

Our PHOTON III family matches these requirements perfectly:

- The largest active area available in an in-house detector with no gaps or dead areas.
- True photon-counting detection for the accurate measurement of weak reflections.
- Operation in charge-integration mode eliminates pulse pileup and charge sharing noise.

The latest generation of detectors for synchrotrons and XFELs uses a new, improved technique for photon counting based on charge integration.


These new charge-integrating detectors, such as the JUNGFRÄU (SWISSFEL) and the ePIX (LCLS), do not suffer anymore from pulse pileup or charge-sharing noise and thus have proven to produce markedly superior data compared to conventional Hybrid Photon-Counting (HPC) detectors.

The PHOTON III is the first in-house detector based on this latest charge-integration technology. The PHOTON III thus features zero detector noise and also eliminates charge sharing and pulse pileup noise that plague conventional HPCs.

The PHOTON III also features an advanced rare-earth X-ray converter that has up to three times the Detective Quantum Efficiency of silicon sensors.

The unique combination of these features makes the PHOTON III the best in-house detector.

The PHOTON III employs large monolithic sensors and therefore has no dead areas or blind gaps: 100% of the surface is fully sensitive to X-rays.



Lysozyme electron density map from the PHOTON III dataset refined to 1.55 Å, contoured at 1σ.

PHOTON III: Large active area – no gaps

- Largest active area for highest redundancy
- No gaps or dead areas
- Charge integration for best linearity
- Photon-counting for best sensitivity
- High dynamic range
- All air-cooled



Best Data

The PHOTON III is the only laboratory detector featuring the latest charge-integrating, photon-counting technology. It excels at both weak and strong reflections, ensuring superior data.



Highest Quality

Three-year warranty, air-cooling and no maintenance go hand in hand to deliver a long-lasting, highly reliable detector.




No blind areas

PHOTON III detectors feature large monolithic sensors. Even multi-modular versions do not suffer from dead areas between the sensors. This allows faster data collection and greater completeness.



Outstanding linearity

Charge integration means no pulse pileup noise and thus the best linearity for strong reflections.



The photon-counting PHOTON III detector delivers outstanding data quality from even the most challenging samples.



Large Active Area

The PHOTON III offers the largest active area for the home laboratory – letting you capture more reflections in just one detector setting and get better separation of closely spaced reflections.



High Speed

With a detector frame rate of 70 Hz, zero readout dead time and shutterless operation, data are acquired quickly and accurately.




Charge Integration

Charge integrating detectors record a very high speed 'movie' of the diffraction pattern. A fast processor then looks for X-ray 'hits' in each frame of the movie. The result: photon counting with no charge sharing or pulse pileup noise for optimal data quality.

Love KAPPA geometry? Perfect, we invented it!

- PHOTON III detector
photon-counting
CPAD technology
- Large, up to 280 cm²
active area, no gaps


PHOTON III
detector




- Motorized detector
distance for best
spatial resolution
- Optimized 2theta
range for highest
diffraction angles

- KAPPA goniometer for
greatest flexibility
- Smallest sphere of
confusion


ISX Stage,
automated in situ
crystal screening



I μ S DIAMOND II
microfocus tube



METALJET D2 PLUS
liquid metal source



- METALJET liquid metal source for highest beam intensity
- Excellent uptime and beam stability

Bruker's four-axis KAPPA goniometer features an ergonomic, open-geometry design, offering:

- Easy sample mounting, harvesting, and monitoring.
- Automatic detector-to-sample distance based on unit cell dimensions and crystal quality.
- Large-unit-cell samples with the shortest detector-to-sample distance.
- Friedel pairs on the same frame.
- Optimum performance for in-house S-SAD phasing experiments.
- Ultimate sample positioning freedom, for unconstrained multiplicity of observations.
- Extremely high angular precision and fast goniometer positioning for productivity.

The same inspiration behind the KAPPA goniometer has also produced a new beam path. Designed for the easiest downstream alignment, and maximizing the accessible 2theta range for high-resolution work, the new beam path allows data resolution to the angstrom – or even better. The beam path design is mirrored by our new real-time path planning software, which combines 3-D models of the current hardware configuration with advanced trajectory algorithms for unprecedented goniometer control. Extensions with custom 3-D models are also possible, which makes it easy to add individual hardware.

The D8 QUEST and D8 VENTURE are both built for excellent sample access and visibility, and Bruker's KAPPA goniometer allows for easy sample mounting and retrieval. The absolutely open design protects your investment with maximum flexibility for future extensions.

D8 Structural Biology Solutions – an abundance of unparalleled innovations.

D8 QUEST and D8 VENTURE: Systems as individual as your research

Sample mounting and alignment are easy and straightforward: large doors give you excellent access to the goniometer, and LED illumination guarantees perfect visibility of your crystal. Polarizing lenses and dimmable LED sample illumination greatly enhance crystal image quality. The crystal image is captured by a high-resolution video microscope and can be viewed remotely throughout the experiment.

With our D8 Macromolecular Solutions, we offer a pioneering diffractometer concept with flexibility and modularity. The D8 QUEST and D8 VENTURE can be perfectly configured for the demands of any imaginable application in protein single crystal X-ray diffraction.

Just follow three simple steps:

1. Choose your preferred X-ray source and optics and combine it with the best KAPPA goniometer.
2. Add the state-of-the-art photon-counting PHOTON III detector with an active area fitting your needs.
3. Complete your system with accessories, such as the ISX stage or the SCOUT sample changing robot.

Now you are ready to run the most sophisticated experiments.

While the experiment is in progress, your sample is constantly exposed to X-rays and continuously rotated, maximizing data acquisition efficiency and optimizing data quality. Finally, the system's firmware provides the real-time status of the configuration and components that is not only used for efficient data collection planning, but is also carefully logged and safely stored with all other experiment information, ready to use for automated report generation and publications from your most sophisticated structure experiments.


D8 QUEST and D8 VENTURE – a pioneering concept of flexibility and modularity.

D8 QUEST

- Small footprint without compromises
- Accommodates μ S DIAMOND II configurations
- Exterior dimensions: 187 cm \times 130 cm \times 114 cm (h \times w \times d)

Configuration example:

- μ S DIAMOND II source
- HELIOS optics
- KAPPA goniometer
- PHOTON III detector





D8 VENTURE

- More room for more experimental flexibility
- Accommodates μ S DIAMOND II or METALJET D2 PLUS configurations
- Exterior dimensions: 202 cm \times 168 cm \times 129 cm (h \times w \times d)

Configuration example:

- METALJET D2 PLUS source
- HELIOS optics
- KAPPA goniometer
- PHOTON III detector



Automation tools: Unleash the full power of your D8 System

Make your sample get to the point!

Small crystals and small X-ray beams require ultra-accurate centering, which pushes mechanical goniometer heads to their limits with respect to accuracy, precision, and backlash. Our automated goniometer head (AGH) is a masterpiece of electronics and mechanical engineering with its high level of miniaturization, precision, and reliability. It allows fully automatic centering of your magnetic-pin-mounted sample.

SCOUT sample changer: high throughput meets high reliability.


The SCOUT sample changing robot combines proven hardware with a powerful user interface. A robust piezo-driven automated goniometer head ensures accurate automated centering of the sample in the X-ray beam. The most important feature of SCOUT is its outstanding reliability. SCOUT meets all the latest CE safety certification requirements for automated operation.

ISX stage: in-situ plate screening made easy.


Many researchers are now interested in screening crystal diffraction quality within the crystallization plate. The ISX stage is a powerful new tool, both for screening crystals and for collecting room-temperature data using the latest serial crystallography techniques – recently pioneered at synchrotron beamlines. The ISX is completely motorized. It mounts easily onto the KAPPA goniometer using a kinematic mount. So it can be installed within minutes and then simply removed after use to continue standard operation with cryo-cooled samples.



AGH for fully automated sample centering




SCOUT: Automated cryo-cooled protein crystal handling



ISX Stage: Automated in-situ crystal screening



PROTEUM



Ultimate ease-of-use, highest flexibility – PROTEUM

Our PROTEUM software guides you through the entire experiment with minimum input and maximum graphical feedback and provides the most complete semi-automated pipeline: suggesting proper defaults whenever possible and asking for your expert decision whenever necessary. From quality assessment to phasing, PROTEUM is easy to use and lets you launch the underlying modules with a single mouse click.

Stay informed about the progress and quality of the experiment via the intuitive GUI. Use the world's best engines – all included in the suite – for crystal screening, data acquisition, data integration, and scaling to generate the best data for phasing.


Benefit from built-in expert knowledge about instrument geometry and data collection strategies.



Crystal Screening

Preset or user-defined crystal mounting positions combined with all required tools for frame analysis convert crystal screening from a tedious chore to a real pleasure.


Get fast and reliable feedback on the diffraction limits via resolution rings and spot separation.



Crystal Snapshots and Movies

The excellent crystal illumination and high-resolution camera enable screen shots from the sample under investigation.


Pictures can be used for internal documentation, teaching or publication... or let the software take a movie.




Unit Cell Determination

An ingenious combination of fast Fourier and difference vector techniques indexes the most difficult data with absolute reliability.

Full nonlinear least-squares cell refinement with graphical feedback, overlay of spots, Bravais lattice determination, and tools for easy matrix manipulation complete this module.






Data Integration

True 3-D integration with algorithms optimized for narrow scans.

Includes the display of integration progress and quality, 3-D reflection profiles, multiple-component spot overlays, and many more.


The best software with on-the-fly integration and extensive feedback.



PROTEUM.SCREENING

An easy-to-use interface for screening experiments with automated scoring feedback drives the SCOUT robot.

The GUI makes it easy to select the best crystal for synchrotron data collection – or let the plug-in collect the data for you.



ISX Stage

The ISX software is optimized for maximum productivity and ease of use. Job queuing allows for unattended screening of multiple samples.

Overview of Features and Benefits

	Feature		Benefit
PHOTON III Detector	Photon-counting pixel-array detector	5th generation pixel array detector technology	Higher speed and sensitivity, best data quality
	Large active area	100 × 70 mm ² 100 × 140 mm ² 200 × 140 mm ²	Higher data redundancy, minimized background scatter
	Fast readout time	14 msec	Faster data collection
	Very high countrate	Up to 4 × 10 ⁶ counts/pixel-sec	Superior data precision
	Readout dead time	0 sec (full shutterless operation)	
	Monolithic silicon sensors		No gaps, no dead areas
	High dynamic range		No reflection overloads
	Mixed mode, integrating detection		More accurate reflection intensities
	Small pixel size, minimized point-spread		Best long axis resolution, best spot separation
	High count-rate linearity	< 1% nonlinearity up to full count rate	Better <i>R</i> -factors
	High detective quantum efficiency (DQE)		
	No operating gas or cooling water	Completely sealed design, air-cooled	No maintenance, high uptime
	High reliability	Warrantied for 3 years	
IμS DIAMOND II Source	Very high intensity beam, completely air-cooled, revolutionary cooling technology	Cu radiation	Rotating performance from a microfocus sealed tube, ultra-low maintenance, optional dual-wavelength configuration
METALJET D2 PLUS Source	Highest intensity beam	Ga radiation	Ultimate performance for small, weakly diffracting samples
D8 Goniometer	Very low sphere of confusion	< 7 μ m	Best data quality
	Kappa geometry	Highest flexibility	Easy sample mounting
	High speed	Up to 1,200 deg/min (omega)	Faster data collection
SCOUT Sample changer	Reliable, small-footprint cryo-cooled protein crystal handling unit for automated crystal screening and data collection	Six-axis robot with auto-refill sample dewar	Enhances system's throughput, ideal for identifying best crystals before a synchrotron trip
ISX Screening stage	Versatile stage for in-situ plate screening	Compatible with all SBS-format multi-well plates, with access to all wells in one setting. Identifies your best crystals, and collects complete data sets at room temperature.	

Bruker AXS
info.baxs@bruker.com

Worldwide offices
bruker.com/baxs-offices

Online information
bruker.com/scd

